

NA-Online-Toolbox: Example Signal-Processing Chain

Datum:	03.03.2015
Autor:	Christof Pieloth
Änderungen:	2015-03-03: First Version

Table of Contents

Remarks.....	2
Offline/Simulated EEG Streaming from FIF File.....	3
Online EEG Streaming with FieldTrip Buffer.....	4
Online MEG Streaming with FieldTrip Buffer (in progress).....	6

Remarks

- Save a signal-processing chain: *File* → *Save Project* → *Save Modules Only*
- Load a signal-processing chain: *File* → *Dataset or Project*
- Confirm changes in input boxes by pressing *Enter* on keyboard, otherwise changes are not applied!
- Some changes are not applied during a running measurement or just from the next block/packet
- FieldTrip Buffers are located in: `fieldtrip/realtime/src/acquisition`
You may have to compile it for your DAQ system.

Offline/Simulated EEG Streaming from FIF File

Instructions:

- Forward solution and source space must be available!
- Setup the NA-Online toolbox (see below)
- Start the streaming by pressing *Send EMM* in FIFF Reader module
- Depending on the setup, you may have to press *Start streaming* in EMM Streaming module

NA-Online Toolbox:

1. FIFF Reader
 - Load FIF file
 - Optional: load additional data
2. FIFF Data Appender
 - Load forward solution/leadfield
 - Load source space
3. EMM Streaming
 - Set block size, e.g. 1000ms
 - Optional: enable *Auto start*
4. FIR Filter
 - Choose filter and windowing type
 - Set sampling and cutoff frequencies
 - Apply settings with *Calculate Filter Coefficients*
5. Epoch Separation
 - Set pre-/post-samples
 - Set trigger code
 - Set event channel (zero-based indexing)
 - Apply settings with *(Re)set*
6. Epoch Averaging
 - Choose average type: total or moving
 - Set *TBase* for baseline correction
 - Apply settings with *(Re)set*
7. Source Reconstruction
 - Set compute modality to EEG
 - Set weighting type: MN (or WMN, in progress)
 - Set SNR

Online EEG Streaming with FieldTrip Buffer

Instructions:

- Do steps from external components (see below)
- Setup the NA-Online toolbox (see below)
- In FieldTrip Real-time Client module press *Connect* and *Start*
- Start your measurement on acquisition system
- NA-Online toolbox starts processing when receiving first data

External Components:

1. Extract fiducial points of subject in AC-PC coordinates (LPA, RPA and Nasion), e.g. from BEM head/skin layer
2. Generate virtual EEG sensors from subject's BEM with `nao_eeg_sensor_generator <args>` (in `na-online_tools`)
3. Generate a high definition forward solution/leadfield matrix with MNE toolbox: `mne_forward_solution <args>`
4. Setup acquisition system
5. Start FieldTrip Buffer on acquisition system, e.g.:
`./neuromag2ft --bufport 4217`

NA-Online Toolbox:

1. FieldTrip Real-time Client
 - Set connection type
 - Set host IP and port
 - Set block size
2. FIFF Data Appender
 - Load BEM layers
 - Load high definition forward solution/leadfield
 - Load source space
3. Alignment
 - Set fiducial points in AC-PC coordinates
4. Leadfield Interpolation
5. FIR Filter
 - Choose filter and windowing type
 - Set sampling and cutoff frequencies
 - Apply settings with *Calculate Filter Coefficients*
6. Epoch Separation
 - Set pre-/post-samples
 - Set trigger code
 - Set event channel (zero-based indexing)
 - Apply setting with *(Re)set*
7. Epoch Averaging
 - Choose average type: total or moving

- Set *TBase* for baseline correction
 - Apply setting with *(Re)set*
8. Source Reconstruction
- Set compute modality to EEG
 - Set weighting type: MN (or WMN, in progress)
 - Set SNR

Online MEG Streaming with FieldTrip Buffer (in progress)

Instructions:

- Do steps from external components (see below)
- Setup the NA-Online toolbox (see below)
- In FieldTrip Real-time Client module press *Connect* and *Start*
- Start your measurement on acquisition system

External Components:

1. Setup acquisition system
2. Do an initial Head Position measurement for a reference position
3. Generate a forward solution/leadfield matrix for the reference position with MNE toolbox:
`mne_forward_solution <args>`
4. Start FieldTrip Buffer on acquisition system, e.g.:
`./neuromag2ft --bufport 4217`

NA-Online Toolbox:

1. FieldTrip Real-time Client
 - Set connection type
 - Set host IP and port
 - Set block size
2. FIFF Data Appender
 - Load BEM layers
 - Load high definition forward solution/leadfield
 - Load source space
3. Head Position Estimation
 - Note: needs unfiltered continuous data stream, frequency based
 - Set HPI coil frequencies
 - Apply settings for *HPI Signal Extraction*
 - Setup *Head Position Estimation*, e.g. epsilon and initial/guessed transformation
4. FIR Filter
 - Choose filter and windowing type
 - Set sampling and cutoff frequencies
 - Apply settings with *Calculate Filter Coefficients*
5. Epoch Separation
 - Set pre-/post-samples
 - Set trigger code
 - Set event channel (zero-based indexing)
 - Apply setting with *(Re)set*
6. Head Position Correction
 - Note: must be placed before averaging, baseline correction, ...
 - Set movement threshold
 - Set sphere radius for generated simple source space
 - Open FIF file with initial Head Position measurement, *Ref. Position*

7. Epoch Averaging
 - Choose average type: total or moving
 - Set *TBase* for baseline correction
 - Apply setting with *(Re)set*
8. Source Reconstruction
 - Set compute modality to MEG
 - Set weighting type: MN (or WMN, in progress)
 - Set SNR