

Clone - C:\Users\rogera\Documents\Digi_ESP\workspace_75\HG1E V3.x\rapIERx

Source:

Destination:

☒ Options

☒ Clone to revision:

☐ Do not update the new working directory

☐ Use pull protocol to copy metadata

☐ Use uncompressed transfer



☐ Include patch queue

☐ Use proxy server

☐ Do not verify host certificate

☐ Remote command:

Hg command:

rapIERx			
Graph	Rev	Branch	Description
	75+	default	★ Working Directory ★
	75	default	default tip Added tag V2.3 for changeset 5ce5d2f7711c